

Geneious workflow for SARS-CoV-2 genome assembly - Download and installation instructions

1 Introduction

This is a workflow for [Geneious Prime](#) that automates the bioinformatics analysis of SARS-CoV-2 genomes from Nanopore data. The workflow runs on FASTQ files and outputs SARS-CoV-2 consensus genome sequences. For instructions on how to perform basecalling on Nanopore data and generate FASTQ files see [Nanopore's instructions](#).

There are two workflows available, “SARS-CoV-2_Assembly_Basic” and “SARS-CoV-2_Assembly_WrapperPlugins”.

- “SARS-CoV-2_Assembly_Basic” performs basic analysis including trimming, mapping, consensus generation and annotation.
- “SARS-CoV-2_Assembly_WrapperPlugins” also includes pango-lineage assignment, using pangolin (O’Toole *et al.*, 2021) and clade assignment, mutations calling and quality checks using Nextclade (Hadfield *et al.*, 2018) as plugins in the workflow.

Both workflows can be used with the ARTIC primers (Quick, 2020) and midnight primers (Freed *et al.*, 2020). The workflows are designed to work on Windows, Linux and Mac OS.

2 Files, Software & Set-up

2.1 Software requirements

Software requirements are presented in Table 1. Some software is required regardless of workflow, while some software is only required for “SARS-CoV-2_Assembly_WrapperPlugins”.

Table 1. A list of software required for a functioning workflow.

Software	Link	Required for
Geneious Prime	geneious.com	*Both workflows
Python	python.org	*Both workflows
Docker Desktop	docker.com	SARS-CoV-2_Assembly_WrapperPlugins
Nextclade CLI	Nextclade-cli	SARS-CoV-2_Assembly_WrapperPlugins

2.2 Setting up the workflows

2.2.1 Download the files from GitHub


Download the GitHub repository by pressing ‘<> Code’ and ‘Download ZIP’

Or use git clone:

```
git clone https://github.com/clinical-genomics-uppsala/Geneious_SARS-CoV-2.git
```

2.2.2 Import data to Geneious

The workflow is designed to work with the structure presented in Figure 1. At the folder overview at the left side in Geneious Prime, select **Local** (or select a subfolder of your choice). Press **Add > Import Folder...** and select the folder “SARS-CoV-2_public_workflows” available in the folder downloaded from GitHub.

A screenshot of the Geneious Prime interface showing a folder tree. The root folder is 'SARS-CoV-2_public_workflows' with 0 items. It contains several subfolders: 'Annotation_features' (0 items), 'Common mutations' (7 items), 'Mutations of interest' (8 items), 'Primer Sars-CoV-2 Artic' (196 items), 'Primer Sars-CoV-2 Midnight1200' (58 items), and 'SARS-CoV-2_reference-sequence' (1 item).

▼	SARS-CoV-2_public_workflows	0
▼	Annotation_features	0
	Common mutations	7
	Mutations of interest	8
	Primer Sars-CoV-2 Artic	196
	Primer Sars-CoV-2 Midnight1200	58
	SARS-CoV-2_reference-sequence	1

Figure 1. An illustration on what the folder structure will look like when imported into Geneious.

2.2.3 Import the workflow to Geneious

Select Tools > *Workflows > Manage Workflows...

Click on **Import** and then select the `.geneiousWorkflow` file that you want to use, “SARS-CoV-2_Assembly_Basic” or “SARS-CoV-2_Assembly_WrapperPlugins”. Once imported, the workflow should appear at the bottom of the list of existing workflows.

Next ensure that the different steps of the workflow are linked to the files imported in section 2.2.2

Double-click on the workflow and the window “Edit Workflow” should appear (Figure 2). The steps that may need manual editing are:

- Trim using BBDuk
- Align/Assemble -> Map to Reference
- Annotate from Database (1)
- Annotate from Database (2)
- Export (last step)




Figure 2. The steps that may require manual modifications have been highlighted. Note that “SARS-CoV-2_Assembly_WrapperPlugins” looks slightly different but contains the same steps that need the same type of modifications.

Trim using BBDuk

Double-click “Trim using BBDuk”. In the window “Edit Trim using BBDuk”, select **Choose...** under **Trim Adapters > Adapters**. Select “Primer Sars-CoV-2 Artic” or “Primer Sars-CoV-2 Midnight1200” depending on which set of primers you plan to use. Select **Apply & Close**.

Options to expose to user when workflow is run

- ☒ Expose no options
- ☐ Expose all options
- ☐ Expose some options

Optionally label exposed options as: ☐ Access exposed options via button

Expose: With Alternative Label:

All Operation Options (those not exposed to workflow user and default values for options that are exposed)

BBDuk Adapter/Quality Trimming Version 38.84 by Brian Bushnell

☒ Trim Adapters

Adapters:

Trim:

Kmer Length:

Maximum Substitutions:

Maximum Substitutions + INDELS:

☐ Trim partial adapters from ends with kmer length:

☒ Trim Low Quality

Trim:

Minimum Quality:

☐ Trim adapters based on paired read overhangs

Minimum Overlap:

☒ Discard Short Reads

Minimum Length: bp

All Operation Options (those not exposed to workflow user and default values for options that are exposed)

☐ Trim Low Complexity (Entropy)

Minimum Entropy:

Entropy Window Size:

Entropy Kmer Size:

☐ Keep original order (slower, but ensures results are reproducible)

Maximum memory to use: MB

Custom BBDuk Options:

Note: The full command line and output from BBDuk is available from the 'Info' tab of the results

*Note: If you plan to use both sets of primers and don't want to manually select the primers set everytime you switch between them, it is possible to create two workflows where each workflow has selected one set of primers each. One way to do this is to go to **Manage Workflows**, select the imported workflow and click **Copy**. You can edit the name of the copied workflow and choose a different set of primers for this workflow.

Align/Assemble -> Map to Reference

The alignment step uses Minimap2. If it is not already installed, go to Tools -> Plugins and under “Available plugins”, the plugin named Minimap should be available for installation. It is also possible to download and install Minimap into Geneious Prime at geneious.com/plugins/minimap2.

Double-click “Align/Assemble -> Map to Reference”. Ensure that the reference sequence has been selected. If it says *automatic*, select **Choose...** and add the reference sequence *NC_045512.fasta*. Ensure that under “Method”, the “Mapper:” has Minimap2 selected. Select **Apply & Close**.

Edit Align/Assemble -> Map to Reference

Align sequences or reads to a reference. Can be used for re-assembly, variant finding, locating a sub-sequence etc

Options to expose to user when workflow is run

- ☒ Expose no options
- ☐ Expose all options
- ☐ Expose some options

Optionally label exposed options as: ☐ Access exposed options via button ?

Expose: With Alternative Label: + -

All Operation Options (those not exposed to workflow user and default values for options that are exposed)

Data

- ☒ Dissolve contigs and re-assemble

Reference Sequence: ?

☐ Assemble by: part of name, separated by

☒ Assemble each sequence list separately

Method

Mapper: ?

[Minimap 2.24](#)

Data Type: ?

☐ Override K-Mer Length: ?

☒ Include Secondary Alignments

Maximum secondary alignments per read: ?

Minimum secondary to primary alignment score ratio:

Note: Paired reads can be set up or changed using **Sequence > Set Paired Reads**

Trim Before Mapping

- ☐ Use existing trim regions
- ☒ Remove existing trim regions from sequences
- ☐ Re-trim sequences
- ☐ Do not trim (discard trim annotations)

Results

Assembly Name ...

- ☐ Save assembly report
- ☐ Save list of unused reads
- ☐ Save list of used reads ☐ Include mates
- ☐ Save in sub-folder
- ☒ Save contigs
- ☐ Save consensus sequences

More Options

Annotate from Database (1)

Double-click the first out of the two steps called “Annotate from Database”. In the window “Edit Annotate from Database”, click on **Reference Features** (located right of **Source:**). Go to the folder *Annotation_Features* and select the subfolder *Common mutations*. Select **Apply & Close**.

Edit Annotate from Database

Annotate sequence with specified annotations based on similarity

Options to expose to user when workflow is run

- ☒ Expose no options
- ☐ Expose all options
- ☐ Expose some options

Optionally label exposed options as: ☐ Access exposed options via button

Expose: Source: With Alternative Label:

All Operation Options (those not exposed to workflow user and default values for options that are exposed)

Source: Common mutations

Similarity: 100 %

Find: ☐ Best match ☒ All matching annotations

All Operation Options (those not exposed to workflow user and default values for options that are exposed)

Source database

Include from source folder:

- ☒ Nucleotide sequences
- ☐ Protein sequences (query translated and compared in all six frames)
- ☒ Sequences in subfolders
- ☐ Truncated annotations
- ☐ Unannotated sequences (transferred as Misc Feature type annotations)

Search parameters

- ☐ Search for Gateway sites
- ☐ Adjust CDS boundaries up to 50 bp to match nearest ORF
- Best Match Criteria:
- Exclude annotations of same type that overlap with best match by 80 %
- Index Length:
- Nucleotide Search: 5
- Translation Search: 4

Types

☐ Only Annotate:

Gene

Coverage - High

Coverage - Low

CRISPR.

Fewer Options ^ Reset to Defaults **Apply & Close** Apply Cancel

Annotate from Database (2)

Double-click the second out of the two steps called “Annotate from Database”. In the window “Edit Annotate from Database”, click on **Reference Features** (located right of **Source:**). Go to the folder *Annotation_Features* and select the subfolder *Mutations of interest*. Select **Apply & Close**.

Edit Annotate from Database

Annotate sequence with specified annotations based on similarity

Options to expose to user when workflow is run

- ☒ Expose no options
- ☐ Expose all options
- ☐ Expose some options

Optionally label exposed options as: ☐ Access exposed options via button

Expose: With Alternative Label:

All Operation Options (those not exposed to workflow user and default values for options that are exposed)

Source:

Similarity: 80 %

Find: ☐ Best match

☒ All matching annotations

All Operation Options (those not exposed to workflow user and default values for options that are exposed)

Source database

Include from source folder:

- ☒ Nucleotide sequences
- ☐ Protein sequences (query translated and compared in all six frames)
- ☒ Sequences in subfolders
- ☒ Truncated annotations
- ☐ Unannotated sequences (transferred as Misc Feature type annotations)

Search parameters

- ☐ Search for Gateway sites
- ☐ Adjust CDS boundaries up to bp to match nearest ORF

Best Match Criteria:

Exclude annotations of same type that overlap with best match by %

Index Length:

Nucleotide Search:

Translation Search:


Types

☐ Only Annotate:

Coverage - High
Coverage - Low
CRISPR.
...

Export (last step)

Double click the step “Export” located at the end of the workflow. Select **Browse** next to the option **Export to Folder**. Select a folder on your computer where you want your final results of the workflow to be located. Select **Apply & Close**.



2.3 Add the plugins for the “SARS-CoV-2_Assembly_WrapperPlugins.geneiousWorkflow”

If you are using the basic workflow without the pango-lineage assignment plugin (*SARS-CoV-2_Assembly_Basic.geneiousWorkflow*), you can skip ahead to section “3 Running the workflow in Geneious”.

If you want to use the workflow called “SARS-CoV-2_Assembly_WrapperPlugins.geneiousWorkflow” some additional steps are necessary.

2.3.1 Import the plugins to Geneious

The plugins with installation instructions are available at these repositories:

- [pangolin](#)
- [Nextclade](#)

Select the relevant plugins depending on your operating system. Ensure that the plugins are successfully imported into Geneious. They should be visible under **Tools**.

2.3.2 Add the plugins to the workflow

Select **Tools > Workflows > Manage Workflows...**

Double-click the workflow “SARS-CoV-2_Assembly_WrapperPlugins.geneiousWorkflow”.

A box with all the steps of the workflow should be displayed.

Select the step that’s below to where the plugin should be (Figure 3) and press **+ Add Step**.

The plugins can be found if selecting **Add Operation (from x available)...**

The two plugin steps should now be located in the correct part of the workflow (Figure 3). If the steps are positioned in the wrong order, you can drag them to the correct position.





Figure 3. A screenshot of the Editing mode of “SARS-CoV-2_Assembly_WrapperPlugins.geneiousWorkflow”. The locations of the plugin steps have been highlighted.

The plugins have some user options that need to be specified. This is done by double clicking the plugin step and a new box will be displayed.

For the pangolin plugin:

Docker path: The path to docker installation source.

Geneious data path: Set this path to the folder *transient* in your local Geneious directory.



Edit Pangolin (linux/mac)

Options to expose to user when workflow is run

- ☒ Expose no options
- ☐ Expose all options
- ☐ Expose some options

Optionally label exposed options as: ☐ Access exposed options via button ?

Expose: options.Docker path With Alternative Label: (Optional) + -

All Operation Options (those not exposed to workflow user and default values for options that are exposed)

Docker path


Geneious data path

More Options ▾ **Apply & Close** Apply Cancel

For the Nextclade plugin:

Nextclade path: The path to the Nextclade installation source.

Path for output files: This is a path of your choice, where you'd like the Nextclade output files to be located.



Edit Nextclade (linux/mac)

Options to expose to user when workflow is run

- ☒ Expose no options
- ☐ Expose all options
- ☐ Expose some options

Optionally label exposed options as: ☐ Access exposed options via button ?

Expose: options.Nextclade path With Alternative Label: (Optional) + -

All Operation Options (those not exposed to workflow user and default values for options that are exposed)

Nextclade path

Path for output files

More Options ▾ **Apply & Close** Apply Cancel

Finally, a folder where to export the result files of the two plugins should be specified. Two positions under the pangolin step, there is a step called “Export”. Double click on it. In the window “Edit Export”, go to the line “Export to Folder” and select **Browse** to specify the output folder of your results. Repeat the same procedure with the step “Export”, located two positions under the Nextclade step.

3 Running the workflow in Geneious

3.1 Run the workflow

Select the FASTQ file(s) to analyse. Then proceed to **Tools > Workflows > “name of the workflow”** and press **OK**.

3.2 Output data

Once the workflow has completed running, new output files should have been generated.

Three subfolders have now been created from where the input FASTQ-files were located:

- Read Mappings
- Consensus
- Reports (only plugin workflow)

The consensus file(s) in fasta-format will also have been exported to the folder selected in the workflow.

For the “SARS-CoV-2_Assembly_WrapperPlugins.geneiousWorkflow”, there should also be:

- X Sequences Pangolin.txt
- X Sequences Nextclade.txt
- a folder called “nextclade_final_<today’s date>” containing Nextclade’s output files according to [Nextclade's output files](#).

References

- Freed, N.E. *et al.* (2020) 'Rapid and inexpensive whole-genome sequencing of SARS-CoV-2 using 1200 bp tiled amplicons and Oxford Nanopore Rapid Barcoding', *Biology Methods & Protocols*, 5(1), p. bpaa014. Available at: <https://doi.org/10.1093/biomethods/bpaa014>.
- Hadfield, J. *et al.* (2018) 'Nextstrain: real-time tracking of pathogen evolution', *Bioinformatics (Oxford, England)*, 34(23), pp. 4121–4123. Available at: <https://doi.org/10.1093/bioinformatics/bty407>.
- O'Toole, Á. *et al.* (2021) 'Assignment of epidemiological lineages in an emerging pandemic using the pangolin tool', *Virus Evolution*, 7(2), p. veab064. Available at: <https://doi.org/10.1093/ve/veab064>.
- Quick, J. (2020) 'nCoV-2019 sequencing protocol v3 (LoCost)'. Available at: <https://www.protocols.io/view/ncov-2019-sequencing-protocol-v3-locost-bh42j8ye> (Accessed: 29 March 2023).